

CD58 Rabbit mAb

Catalog No: #56951

Package Size: #56951-1 50ul #56951-2 100ul

Description

Product Name	CD58 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP FC
Species Reactivity	Human
Specificity	CD58 Antibody detects endogenous levels of total CD58
Immunogen Description	A synthesized peptide derived from human CD58
Conjugates	Unconjugated
Other Names	AG3; CD58; LFA3;
Accession No.	Uniprot:P19256
Calculated MW	50-75kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IP:1:50FC:1:50

Note: This product is for in vitro research use only and is not intended for use in humans or animals.