

ARPC2 Rabbit mAb

Catalog No: #56974

Package Size: #56974-1 50ul #56974-2 100ul


Description

Product Name	ARPC2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	ARPC2 Antibody detects endogenous levels of total ARPC2
Immunogen Description	A synthesized peptide derived from human ARPC2
Conjugates	Unconjugated
Other Names	ARC34; Arpc2; p34Arc; PNAS139; PRO2446;
Accession No.	Uniprot:O15144
Calculated MW	34kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50FC:1:50

Images



Western blot analysis of ARPC2 expression in 293 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.