

MYL12B Rabbit mAb

Catalog No: #57015

Package Size: #57015-1 50ul #57015-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MYL12B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	MYL12B Antibody detects endogenous levels of total MYL12B
Immunogen Description	A synthesized peptide derived from human MYL12B
Conjugates	Unconjugated
Other Names	MLC 2A; MLC20; MRLC2; MYL12B; MYLC2B;
Accession No.	Uniprot:O14950
Calculated MW	18kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200

Images



Western blot analysis of MYL12B expression in Jurkat cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.