

CLIP170 Rabbit mAb

Catalog No: #57074

Package Size: #57074-1 50ul #57074-2 100ul

Description

Product Name	CLIP170 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	CLIP170 Antibody detects endogenous levels of total CLIP170
Immunogen Description	A synthesized peptide derived from human CLIP170
Conjugates	Unconjugated
Other Names	CLIP; CYLN1; Restin; RSN;
Accession No.	Uniprot:P30622
Calculated MW	130,170kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.