

Syndecan 1 Rabbit mAb

Catalog No: #59388

Package Size: #59388-1 50ul #59388-2 100ul


Description

Product Name	Syndecan 1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	Syndecan 1 Antibody detects endogenous levels of total Syndecan 1
Immunogen Description	A synthesized peptide derived from human Syndecan 1
Other Names	CD138; CD138 antigen; heparan sulfate proteoglycan fibroblast growth factor receptor; Sdc1; SYND1; Syndecan 1; syndecan proteoglycan 1;
Accession No.	Uniprot:P18827
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:100~1:500 ICC/IF 1:50~1:200

Images



Western blot analysis of Syndecan 1 expression in Ramos cell lysate.

Background

Cell surface proteoglycan that bears both heparan sulfate and chondroitin sulfate and that links the cytoskeleton to the interstitial matrix.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.