

ADFP Rabbit mAb

Catalog No: #59404

Package Size: #59404-1 50ul #59404-2 100ul


Description

Product Name	ADFP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	ADFP Antibody detects endogenous levels of total ADFP
Immunogen Description	A synthesized peptide derived from human ADFP
Other Names	ADFP; Adipophilin; Adipose differentiation related protein; ADRP; Perilipin-2; PLIN2;
Accession No.	Uniprot:Q99541
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of ADFP expression in HepG2 cell lysate.

Background

May be involved in development and maintenance of adipose tissue.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.