

ERp57 Rabbit mAb

Catalog No: #59410

Package Size: #59410-1 50ul #59410-2 100ul


Description

Product Name	ERp57 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human
Specificity	ERp57 Antibody detects endogenous levels of total ERp57
Immunogen Description	A synthesized peptide derived from human ERp57
Other Names	ERp57; ERp60; ERp61; Glucose Regulated Protein 58 Kd; HsT17083; p58; PDIA3; Phospholipase C alpha; PI PLC;
Accession No.	Uniprot:P30101
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 IP 1:50

Images



Western blot analysis of ERp57 expression in HepG2 cell lysate.

Background

Detected in the flagellum and head region of spermatozoa (at protein level).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.