

NFAT2 Rabbit mAb

Catalog No: #59437

Package Size: #59437-1 50ul #59437-2 100ul


Description

Product Name	NFAT2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human
Specificity	NFAT2 Antibody detects endogenous levels of total NFAT2
Immunogen Description	A synthesized peptide derived from human NFAT2
Other Names	cytoplasmic 1; NF ATc; NF ATc1; NF-ATc; NF-ATc1; NFAC1; NFAT 2; NFAT2; NFATC 1; NFATc; NFATc1;
Accession No.	Uniprot:O95644
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of NFAT2 expression in Ramos cell lysate.

Background

NFAT proteins are predominantly expressed in cells of the immune system, but are also expressed in skeletal muscle, keratinocytes, and adipocytes, regulating cell differentiation programs in these cells.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.