

## TGM2 Rabbit mAb

Catalog No: #59449

Package Size: #59449-1 50ul #59449-2 100ul


## Description

|                       |                                                                                                                                  |
|-----------------------|----------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | TGM2 Rabbit mAb                                                                                                                  |
| Host Species          | Rabbit                                                                                                                           |
| Clonality             | Monoclonal                                                                                                                       |
| Isotype               | Rabbit IgG                                                                                                                       |
| Purification          | Affinity-chromatography                                                                                                          |
| Applications          | WB IHC IP                                                                                                                        |
| Species Reactivity    | Human Mouse                                                                                                                      |
| Specificity           | TGM2 Antibody detects endogenous levels of total TGM2                                                                            |
| Immunogen Description | A synthesized peptide derived from human TGM2                                                                                    |
| Other Names           | TGM2; Tissue transglutaminase; Transglutaminase C; TG(C); TGC; TGase C; ransglutaminase H; TGase H; Transglutaminase-2; TGase-2; |
| Accession No.         | Uniprot:P21980                                                                                                                   |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.                               |
| Storage               | Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.                                                   |

## Application Details

WB 1:1000~1:5000 IHC 1:50~1:200 IP 1:50

## Images



Western blot analysis of TGM2 expression in HUVEC cell lysate.

## Background

Catalyzes the cross-linking of proteins and the conjugation of polyamines to proteins.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.