

PCSK9 Rabbit mAb

Catalog No: #59450

Package Size: #59450-1 50ul #59450-2 100ul


Description

Product Name	PCSK9 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human
Specificity	PCSK9 Antibody detects endogenous levels of total PCSK9
Immunogen Description	A synthesized peptide derived from human PCSK9
Other Names	HCHOLA3; LDLQCQ1; NARC1; Neural apoptosis regulated convertase 1; PC9; PCSK9; Proprotein convertase 9;
Accession No.	Uniprot:Q8NBP7
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IP 1:20

Images



Western blot analysis of PCSK9 expression in HepG2 cell lysate.

Background

May be implicated in the differentiation of cortical neurons and may play a role in cholesterol homeostasis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.