

gamma Sarcoglycan Rabbit mAb

Catalog No: #59456

Package Size: #59456-1 50ul #59456-2 100ul


Description

Product Name	gamma Sarcoglycan Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human
Specificity	gamma Sarcoglycan Antibody detects endogenous levels of total gamma Sarcoglycan
Immunogen Description	A synthesized peptide derived from human gamma Sarcoglycan
Other Names	35DAG; DAGA4; DMDA; Gamma-sarcoglycan; LGMD2C; MAM; Sarcoglycan gamma; SCARMD2; SCG3; SGCG; TYPE;
Accession No.	Uniprot:Q13326
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 IP 1:50

Images



Western blot analysis of gamma Sarcoglycan expression in Human skeletal muscle lysate.

Background

Component of the sarcoglycan complex, a subcomplex of the dystrophin-glycoprotein complex which forms a link between the F-actin cytoskeleton and the extracellular matrix.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.