

delta Sarcoglycan Rabbit mAb

Catalog No: #59459

Package Size: #59459-1 50ul #59459-2 100ul


Description

Product Name	delta Sarcoglycan Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human
Specificity	delta Sarcoglycan Antibody detects endogenous levels of total delta Sarcoglycan
Immunogen Description	A synthesized peptide derived from human delta Sarcoglycan
Other Names	35DAG; CMD1L; DAGD; Delta-sarcoglycan; LGMD2F; SG delta; SGCD; SGCDP;
Accession No.	Uniprot:Q92629
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 IP 1:50

Images



Western blot analysis of delta Sarcoglycan expression in Human fetal heart lysate.

Background

Component of the sarcoglycan complex, a subcomplex of the dystrophin-glycoprotein complex which forms a link between the F-actin cytoskeleton and the extracellular matrix.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.