

CRYAB Rabbit mAb

Catalog No: #59463

Package Size: #59463-1 50ul #59463-2 100ul


Description

Product Name	CRYAB Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Rat
Specificity	CRYAB Antibody detects endogenous levels of total CRYAB
Immunogen Description	A synthesized peptide derived from human CRYAB
Other Names	AACRYA; Alpha B crystallin; Alpha-crystallin B chain; CRYA2; Cryab; Crystallin alpha B; CTPP2; Heat shock 20 kD like protein; HspB5;
Accession No.	Uniprot:P02511
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:80

Images



Western blot analysis of CRYAB expression in human fetal heart lysate.

Background

May contribute to the transparency and refractive index of the lens.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.