

BNIP3L Rabbit mAb

Catalog No: #59467

Package Size: #59467-1 50ul #59467-2 100ul


Description

Product Name	BNIP3L Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	BNIP3L Antibody detects endogenous levels of total BNIP3L
Immunogen Description	A synthesized peptide derived from human BNIP3L
Other Names	Adenovirus E1B19k binding protein B5; BNIP3a; BNIP3H; BNIP3L; BNIP3L protein; NIP3 like protein X; NIP3L; NIX;
Accession No.	Uniprot:O60238
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Images



Western blot analysis of BNIP3L expression in K562 cell lysate.

Background

Induces apoptosis. Interacts with viral and cellular anti-apoptosis proteins. Can overcome the suppressors BCL-2 and BCL-XL, although high levels of BCL-XL expression will inhibit apoptosis. May function as a tumor suppressor. Inhibits apoptosis induced by BNIP3.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.