

## BIRC2 Rabbit mAb

Catalog No: #59473

Package Size: #59473-1 50ul #59473-2 100ul


## Description

|                       |                                                                                                         |
|-----------------------|---------------------------------------------------------------------------------------------------------|
| Product Name          | BIRC2 Rabbit mAb                                                                                        |
| Host Species          | Rabbit                                                                                                  |
| Clonality             | Monoclonal                                                                                              |
| Isotype               | Rabbit IgG                                                                                              |
| Purification          | Affinity-chromatography                                                                                 |
| Applications          | WB IHC ICC                                                                                              |
| Species Reactivity    | Human                                                                                                   |
| Specificity           | BIRC2 Antibody detects endogenous levels of total BIRC2                                                 |
| Immunogen Description | A synthesized peptide derived from human BIRC2                                                          |
| Other Names           | C-IAP1; IAP homolog B; Inhibitor of apoptosis protein 2; IAP-2; hIAP-2; hIAP2; API1; IAP2; MIHB; RNF48; |
| Accession No.         | Uniprot:Q13490                                                                                          |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.      |
| Storage               | Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.                          |

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200

## Images



Western blot analysis of BIRC2 expression in HepG2 cell lysate.

## Background

Apoptotic suppressor. The BIR motifs region interacts with TNF receptor associated factors 1 and 2 (TRAF1 and TRAF2) to form an heteromeric complex, which is then recruited to the tumor necrosis factor receptor 2 (TNFR2).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.