

Follistatin Rabbit mAb

Catalog No: #59538

Package Size: #59538-1 50ul #59538-2 100ul


Description

Product Name	Follistatin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Rat
Specificity	Follistatin Antibody detects endogenous levels of total Follistatin
Immunogen Description	A synthesized peptide derived from human Follistatin
Other Names	Activin binding protein; Follistatin; FS; fst; Human follistatin gene;
Accession No.	Uniprot:P19883
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of Follistatin expression in HepG2 cell lysate.

Background

Binds directly to activin and functions as an activin antagonist. Specific inhibitor of the biosynthesis and secretion of pituitary follicle stimulating hormone (FSH).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.