

DNA Polymerase beta Rabbit mAb

Catalog No: #59547

Package Size: #59547-1 50ul #59547-2 100ul



Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	DNA Polymerase beta Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	DNA Polymerase beta Antibody detects endogenous levels of total DNA Polymerase beta
Immunogen Description	A synthesized peptide derived from human DNA Polymerase beta
Other Names	DNA directed DNA polymerase beta; DNA pol beta; DNA polymerase beta; Pol B; Pol beta; POLB; Polymerase (DNA directed) beta;
Accession No.	Uniprot:P06746
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50

Images



Western blot analysis of DNA Polymerase beta expression in A431 cell lysate.

Background

Repair polymerase. Conducts gap-filling DNA synthesis in a stepwise distributive fashion rather than in a processive fashion as for other DNA polymerases. Has a 5'-deoxyribose-5-phosphate lyase (dRP lyase) activity.,P06746,Null,<http://www.uniprot.org/uniprot/P06746>,Null 324215,324215,59548,Null,MUC4 Rabbit mAb,Rabbit,Monoclonal,Null,Rabbit IgG,Affinity-chromatography,WB,WB 1:500~1:2000,Human,MUC4 Antibody detects endogenous levels of total MUC4,,,,,Ascites sialoglycoprotein 1; Ascites sialoglycoprotein 2; Ascites sialoglycoprotein; ASGP; ASGP-1; ASGP-2; Muc4; Mucin 4 cell surface associated; Mucin 4 tracheobronchial; Mucin-4 beta chain; Testis mucin;,Uniprot:Q99102,,Null,Null,120kDa,Null,Null,Null,Null,Null,Null,A synthesized peptide derived from human

MUC4,Null,Null,Null,Null,Null,Null,Null,Null,Null,Null,Null,Null,"Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.