

# PMP22 Rabbit mAb

Catalog No: #59556

Package Size: #59556-1 50ul #59556-2 100ul


## Description

|                       |                                                                                                                        |
|-----------------------|------------------------------------------------------------------------------------------------------------------------|
| Product Name          | PMP22 Rabbit mAb                                                                                                       |
| Host Species          | Rabbit                                                                                                                 |
| Clonality             | Monoclonal                                                                                                             |
| Isotype               | Rabbit IgG                                                                                                             |
| Purification          | Affinity-chromatography                                                                                                |
| Applications          | WB                                                                                                                     |
| Species Reactivity    | Human Rat                                                                                                              |
| Specificity           | PMP22 Antibody detects endogenous levels of total PMP22                                                                |
| Immunogen Description | A synthesized peptide derived from human PMP22                                                                         |
| Other Names           | CMT1A; CMT1E; DSS; GAS3; Growth Arrest Specific 3; HMSNIA; HNPP; Peripheral myelin protein 22; PMP22; Sp110; Trembler; |
| Accession No.         | Uniprot:Q01453                                                                                                         |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.                     |
| Storage               | Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.                                         |

## Application Details

WB 1:500~1:2000

## Images



Western blot analysis of PMP22 expression in 293T cell lysate.

## Background

Involved in growth regulation, and in myelination in the peripheral nervous system.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.