

Myelin oligodendrocyte glycoprotein Rabbit mAb

Catalog No: #59557

Package Size: #59557-1 50ul #59557-2 100ul


Description

Product Name	Myelin oligodendrocyte glycoprotein Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	Myelin oligodendrocyte glycoprotein Antibody detects endogenous levels of total Myelin oligodendrocyte glycoprotein
Immunogen Description	A synthesized peptide derived from human Myelin oligodendrocyte glycoprotein
Other Names	BTN6; BTNL11; MOG alpha 5; MOG alpha 6; MOG AluA; MOG AluB; MOG; MOG Ig AluB; MOGIG2; Myelin oligodendrocyte glycoprotein; NRCLP7;
Accession No.	Uniprot:Q16653
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:100~1:500 ICC/IF 1:100~1:500

Images



Western blot analysis of Myelin oligodendrocyte glycoprotein expression in mouse brain lysate.

Background

Mediates homophilic cell-cell adhesion (By similarity). Minor component of the myelin sheath. May be involved in completion and/or maintenance of the myelin sheath and in cell-cell communication.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.