

USP11 Rabbit mAb

Catalog No: #59585

Package Size: #59585-1 50ul #59585-2 100ul


Description

Product Name	USP11 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	USP11 Antibody detects endogenous levels of total USP11
Immunogen Description	A synthesized peptide derived from human USP11
Other Names	UHX1; USP11;
Accession No.	Uniprot:P51784
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 IP 1:100 FC 1:200

Images



Western blot analysis of USP11 expression in (1) HEK293 cell lysate; (2) Mouse testis lysate.

Background

Protease that can remove conjugated ubiquitin from target proteins and polyubiquitin chains. Inhibits the degradation of target proteins by the proteasome. Plays a role in the regulation of pathways leading to NF-kappa-B activation. Plays a role in the regulation of DNA repair after double-stranded DNA breaks.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.