

SKP2 Rabbit mAb

Catalog No: #59586

Package Size: #59586-1 50ul #59586-2 100ul


Description

Product Name	SKP2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	SKP2 Antibody detects endogenous levels of total SKP2
Immunogen Description	A synthesized peptide derived from human SKP2
Other Names	FBL1; FLB1; FBXL1; MGC1366; SKP2;
Accession No.	Uniprot:Q13309
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of SKP2 expression in (1) Jurkat cell lysate; (2) NIH/3T3 cell lysate.

Background

Substrate recognition component of a SCF (SKP1-CUL1-F-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins involved in cell cycle progression, signal transduction and transcription.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.