

FADD Rabbit mAb

Catalog No: #59590

Package Size: #59590-1 50ul #59590-2 100ul


Description

Product Name	FADD Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC FC
Species Reactivity	Mouse
Specificity	FADD Antibody detects endogenous levels of total FADD
Immunogen Description	A synthesized peptide derived from human FADD
Other Names	FADD protein; GIG3; Growth inhibiting gene 3 protein; Mediator of receptor induced toxicity; MORT1;
Accession No.	Uniprot:Q61160
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 FC 1:50

Images



Western blot analysis of FADD expression in Raw264.7 cell lysate.

Background

Apoptotic adaptor molecule that recruits caspase-8 or caspase-10 to the activated Fas (CD95) or TNFR-1 receptors. The resulting aggregate called the death-inducing signaling complex (DISC) performs caspase-8 proteolytic activation. Active caspase-8 initiates the subsequent cascade of caspases mediating apoptosis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.