

RRM2 Rabbit mAb

Catalog No: #59643

Package Size: #59643-1 50ul #59643-2 100ul


Description

Product Name	RRM2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP
Species Reactivity	Human
Specificity	RRM2 Antibody detects endogenous levels of total RRM2
Immunogen Description	A synthesized peptide derived from human RRM2
Other Names	R2;RR2M;
Accession No.	Uniprot:P31350
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 IP 1:50

Images



Western blot analysis of RRM2 expression in HeLa cell lysate.

Background

Provides the precursors necessary for DNA synthesis. Catalyzes the biosynthesis of deoxyribonucleotides from the corresponding ribonucleotides. Inhibits Wnt signaling.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.