

EB3 Rabbit mAb

Catalog No: #59682

Package Size: #59682-1 50ul #59682-2 100ul


Description

Product Name	EB3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	EB3 Antibody detects endogenous levels of total EB3
Immunogen Description	A synthesized peptide derived from human EB3
Other Names	EB3; EBF3; MAPRE3; RP/EB Family; RP3;
Accession No.	Uniprot:Q9UPY8
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:5000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of EB3 expression in Rat muscle cell lysate.

Background

May be involved in microtubule polymerization, and spindle function by stabilizing microtubules and anchoring them at centrosomes. May play a role in cell migration.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.