

Uteroglobin Rabbit mAb

Catalog No: #59689

Package Size: #59689-1 50ul #59689-2 100ul


Description

Product Name	Uteroglobin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse Rat
Specificity	Uteroglobin Antibody detects endogenous levels of total Uteroglobin
Immunogen Description	A synthesized peptide derived from human Uteroglobin
Other Names	Blastokinin; CC10; CC16; CCPBP; CCSP; SCGB1A; UGB; UP1; Uteroglobin;
Accession No.	Uniprot:P11684
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 FC 1:50

Images



Western blot analysis of Uteroglobin expression in A549 cell lysate.

Background

Binds phosphatidylcholine, phosphatidylinositol, polychlorinated biphenyls (PCB) and weakly progesterone, potent inhibitor of phospholipase A2.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.