

GABA A Receptor alpha 5 Rabbit mAb

Catalog No: #59695

Package Size: #59695-1 50ul #59695-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GABA A Receptor alpha 5 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	GABA A Receptor alpha 5 Antibody detects endogenous levels of total GABA A Receptor alpha 5
Immunogen Description	A synthesized peptide derived from human GABA A Receptor alpha 5
Other Names	GAA 5; GAA5; GABRA 5; Gabra5;
Accession No.	Uniprot:P31644
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200

Images



Western blot analysis of GABA A Receptor alpha 5 expression in 293T cell lysate.

Background

GABA, the major inhibitory neurotransmitter in the vertebrate brain, mediates neuronal inhibition by binding to the GABA/benzodiazepine receptor and opening an integral chloride channel.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.