

SMARCC1/BAF155 Rabbit mAb

Catalog No: #59699

Package Size: #59699-1 50ul #59699-2 100ul


Description

Product Name	SMARCC1/BAF155 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human Rat
Specificity	SMARCC1/BAF155 Antibody detects endogenous levels of total SMARCC1/BAF155
Immunogen Description	A synthesized peptide derived from human SMARCC1/BAF155
Other Names	BAF155; CRACC1; Rsc8; SMARCC1; SRG3; SWI3;
Accession No.	Uniprot:Q92922
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 IP 1:50 FC 1:50

Images



Western blot analysis of SMARCC1 expression in HeLa cell lysate.

Background

Involved in transcriptional activation and repression of select genes by chromatin remodeling (alteration of DNA-nucleosome topology). May stimulate the ATPase activity of the catalytic subunit of the complex.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.