

Wnt2b Rabbit mAb

Catalog No: #59715

Package Size: #59715-1 50ul #59715-2 100ul


Description

Product Name	Wnt2b Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	Wnt2b Antibody detects endogenous levels of total Wnt2b
Immunogen Description	A synthesized peptide derived from human Wnt2b
Other Names	WNT13; Wnt2b; XWNT2;
Accession No.	Uniprot:Q93097
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of Wnt2b expression in HepG2 cell lysate.

Background

Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. May be a signaling molecule which affects the development of discrete regions of tissues. Is likely to signal over only few cell diameters.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.