

# RNA Helicase A Rabbit mAb

Catalog No: #59719

Package Size: #59719-1 50ul #59719-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                                                    |
|-----------------------|----------------------------------------------------------------------------------------------------|
| Product Name          | RNA Helicase A Rabbit mAb                                                                          |
| Host Species          | Rabbit                                                                                             |
| Clonality             | Monoclonal                                                                                         |
| Isotype               | Rabbit IgG                                                                                         |
| Purification          | Affinity-chromatography                                                                            |
| Applications          | WB IHC ICC/IF FC                                                                                   |
| Species Reactivity    | Human Mouse Rat                                                                                    |
| Specificity           | RNA Helicase A Antibody detects endogenous levels of total RNA Helicase A                          |
| Immunogen Description | A synthesized peptide derived from human RNA Helicase A                                            |
| Other Names           | DDX 9; DDX9; dhx9; Leukophysin; LKP; NDH2; NDHII; RHA;                                             |
| Accession No.         | Uniprot:Q08211                                                                                     |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage               | Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.                     |

## Application Details

WB 1:1000~1:5000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:20

## Images



Western blot analysis of RNA Helicase A expression in Jurkat cell lysate.

## Background

Unwinds double-stranded DNA and RNA in a 3' to 5' direction. Alteration of secondary structure may subsequently influence interactions with proteins or other nucleic acids. Functions as a transcriptional activator. Component of the CRD-mediated complex that promotes MYC mRNA stability.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.