

USP9x Rabbit mAb

Catalog No: #59724

Package Size: #59724-1 50ul #59724-2 100ul


Description

Product Name	USP9x Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	USP9x Antibody detects endogenous levels of total USP9x
Immunogen Description	A synthesized peptide derived from human USP9x
Other Names	DFFRX; FAF; FafI; Fam; hFAM; MRX99; Usp9x; X chromosome; X-linked;
Accession No.	Uniprot:Q93008
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of USP9x expression in HeLa cell lysate.

Background

Deubiquitinase involved both in the processing of ubiquitin precursors and of ubiquitinated proteins. May therefore play an important role regulatory role at the level of protein turnover by preventing degradation of proteins through the removal of conjugated ubiquitin.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.