

KAT8 / MYST1 / MOF Rabbit mAb

Catalog No: #59754

Package Size: #59754-1 50ul #59754-2 100ul


Description

Product Name	KAT8 / MYST1 / MOF Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	KAT8 / MYST1 / MOF Antibody detects endogenous levels of total KAT8 / MYST1 / MOF
Immunogen Description	A synthesized peptide derived from human KAT8 / MYST1 / MOF
Other Names	hMOF; MOF; MOZ; MYST1; SAS2; SAS3; YBF2; YBF2/SAS3; ZC2HC8;
Accession No.	Uniprot:Q9H7Z6
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:1000 IHC 1:50~1:200 ICC/IF 1:100~1:500 IP 1:100 FC 1:100

Images



Western blot analysis of KAT8 / MYST1 / MOF expression in HeLa cell lysate.

Background

Histone acetyltransferase which may be involved in transcriptional activation. May influence the function of ATM.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.