

SULT2A1 Rabbit mAb

Catalog No: #59757

Package Size: #59757-1 50ul #59757-2 100ul


Description

Product Name	SULT2A1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	SULT2A1 Antibody detects endogenous levels of total SULT2A1
Immunogen Description	A synthesized peptide derived from human SULT2A1
Other Names	HST; ST2; STD; hSTa; DHEAS; ST2A1; ST2A3; DHEA-ST;
Accession No.	Uniprot:Q06520
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:100~1:500 ICC/IF 1:50~1:200 FC 1:80

Images



Western blot analysis of SULT2A1 expression in HepG2 cell lysate.

Background

Catalyzes the sulfation of steroids and bile acids in the liver and adrenal glands.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.