

Apolipoprotein A II Rabbit mAb

Catalog No: #59780

Package Size: #59780-1 50ul #59780-2 100ul


Description

Product Name	Apolipoprotein A II Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human Rat
Specificity	Apolipoprotein A II Antibody detects endogenous levels of total Apolipoprotein A II
Immunogen Description	A synthesized peptide derived from human Apolipoprotein A II
Other Names	APO A2; Apo AII; APOA 2; ApoA II; APOA2; APOAII; Apolipoprotein A II;
Accession No.	Uniprot:P02652
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of Apolipoprotein A II expression in Human plasma lysate.

Background

May stabilize HDL (high density lipoprotein) structure by its association with lipids, and affect the HDL metabolism.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.