

Biglycan Rabbit mAb

Catalog No: #59781

Package Size: #59781-1 50ul #59781-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Biglycan Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP FC
Species Reactivity	Human Mouse Rat
Specificity	Biglycan Antibody detects endogenous levels of total Biglycan
Immunogen Description	A synthesized peptide derived from human Biglycan
Other Names	BGN; Biglycan; DSPG1; PG S1; PGI; SEMDX; SLRR1A;
Accession No.	Uniprot:P21810
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IP 1:50 FC 1:100

Images



Western blot analysis of Biglycan expression in C2C12 cell lysate.

Background

May be involved in collagen fiber assembly.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.