

Fibrinogen gamma chain Rabbit mAb

Catalog No: #59786

Package Size: #59786-1 50ul #59786-2 100ul


Description

Product Name	Fibrinogen gamma chain Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human
Specificity	Fibrinogen gamma chain Antibody detects endogenous levels of total Fibrinogen gamma chain
Immunogen Description	A synthesized peptide derived from human Fibrinogen gamma chain
Other Names	FGG; Fibrinogen gamma chain;
Accession No.	Uniprot:P02679
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 FC 1:50

Images



Western blot analysis of Fibrinogen gamma chain expression in Human plasma lysate.

Background

Fibrinogen has a double function: yielding monomers that polymerize into fibrin and acting as a cofactor in platelet aggregation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.