

Cullin4A Rabbit mAb

Catalog No: #59792

Package Size: #59792-1 50ul #59792-2 100ul


Description

Product Name	Cullin4A Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human
Specificity	Cullin4A Antibody detects endogenous levels of total Cullin4A
Immunogen Description	A synthesized peptide derived from human Cullin4A
Other Names	Cul4a; Cullin-4A;
Accession No.	Uniprot:Q13619
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:5000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50 FC 1:50

Images



Western blot analysis of Cullin 4a expression in HeLa cell lysate.

Background

Core component of multiple cullin-RING-based E3 ubiquitin-protein ligase complexes which mediate the ubiquitination and subsequent proteasomal degradation of target proteins. As a scaffold protein may contribute to catalysis through positioning of the substrate and the ubiquitin-conjugating enzyme.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.