

Plexin A1 Rabbit mAb

Catalog No: #59797

Package Size: #59797-1 50ul #59797-2 100ul


Description

Product Name	Plexin A1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	Plexin A1 Antibody detects endogenous levels of total Plexin A1
Immunogen Description	A synthesized peptide derived from human Plexin A1
Other Names	NOV; NOV; Plexin-A1; Plexin1; PlexinA1; PLXN1; Plxna1;
Accession No.	Uniprot:Q9UIW2
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of Plexin A1 expression in HUVEC cell lysate.

Background

Plays a role in axon guidance, invasive growth and cell migration.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.