

Villin Rabbit mAb

Catalog No: #59798

Package Size: #59798-1 50ul #59798-2 100ul


Description

Product Name	Villin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse Rat
Specificity	Villin Antibody detects endogenous levels of total Villin
Immunogen Description	A synthesized peptide derived from human Villin
Other Names	VIL; VIL1; Villin1;
Accession No.	Uniprot:P09327
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:5000 FC 1:50

Images



Western blot analysis of Villin expression in Caco-2 cell lysate.

Background

Ca(2+)-regulated actin-binding protein.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.