

RhoGDI Rabbit mAb

Catalog No: #59802

Package Size: #59802-1 50ul #59802-2 100ul


Description

Product Name	RhoGDI Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	RhoGDI Antibody detects endogenous levels of total RhoGDI
Immunogen Description	A synthesized peptide derived from human RhoGDI
Other Names	ARHGDI; GDIA1; RHOGDI; RHOGDI-1;
Accession No.	Uniprot:P52565
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200

Images



Western blot analysis of RhoGDI expression in HeLa cell lysate.

Background

Regulates the GDP/GTP exchange reaction of the Rho proteins by inhibiting the dissociation of GDP from them, and the subsequent binding of GTP to them.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.