

Pumilio 1 Rabbit mAb

Catalog No: #59803

Package Size: #59803-1 50ul #59803-2 100ul


Description

Product Name	Pumilio 1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	Pumilio 1 Antibody detects endogenous levels of total Pumilio 1
Immunogen Description	A synthesized peptide derived from human Pumilio 1
Other Names	HsPUM; PUM L1; Pum1; PUMH 1; PUMH; PUMH1; Pumilio1; PUML1;
Accession No.	Uniprot:Q14671
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50 FC 1:60

Images



Western blot analysis of Pumilio 1 expression in HeLa cell lysate.

Background

Mediates post-transcriptional repression of transcripts via different mechanisms: acts via direct recruitment of the CCR4-POP2-NOT deadenylase leading to translational inhibition and mRNA degradation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.