

SENP1 Rabbit mAb

Catalog No: #59805

Package Size: #59805-1 50ul #59805-2 100ul


Description

Product Name	SENP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human
Specificity	SENP1 Antibody detects endogenous levels of total SENP1
Immunogen Description	A synthesized peptide derived from human SENP1
Other Names	SENP1; SuPr2;
Accession No.	Uniprot:Q9P0U3
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:30

Images



Western blot analysis of SENP1 expression in U87-MG cell lysate.

Background

Protease that catalyzes two essential functions in the SUMO pathway: processing of full-length SUMO1, SUMO2 and SUMO3 to their mature forms and deconjugation of SUMO1, SUMO2 and SUMO3 from targeted proteins.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.