

Nucleostemin Rabbit mAb

Catalog No: #59807

Package Size: #59807-1 50ul #59807-2 100ul


Description

Product Name	Nucleostemin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human
Specificity	Nucleostemin Antibody detects endogenous levels of total Nucleostemin
Immunogen Description	A synthesized peptide derived from human Nucleostemin
Other Names	E2IG3; GNL3; NNP47; NS; Nucleostemin;
Accession No.	Uniprot:Q9BVP2
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of Nucleostemin expression in HeLa cell lysate.

Background

May be required to maintain the proliferative capacity of stem cells and may play an important role in tumorigenesis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.