

PICK1 Rabbit mAb

Catalog No: #59812

Package Size: #59812-1 50ul #59812-2 100ul


Description

Product Name	PICK1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP FC
Species Reactivity	Human Mouse Rat
Specificity	PICK1 Antibody detects endogenous levels of total PICK1
Immunogen Description	A synthesized peptide derived from human PICK1
Other Names	PICK 1; Pick1; PRKCABP;
Accession No.	Uniprot:Q9NRD5
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IP 1:50 FC 1:50

Images



Western blot analysis of PICK1 expression in HeLa cell lysate.

Background

Probable adapter protein that bind to and organize the subcellular localization of a variety of membrane proteins containing some PDZ recognition sequence. Involved in the clustering of various receptors, possibly by acting at the receptor internalization level. Plays a role in synaptic plasticity by regulating the trafficking and internalization of AMPA receptors.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.