

p60 CAF1 Rabbit mAb

Catalog No: #59846

Package Size: #59846-1 50ul #59846-2 100ul

Description


Product Name	p60 CAF1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	p60 CAF1 Antibody detects endogenous levels of total p60 CAF1
Immunogen Description	A synthesized peptide derived from human p60 CAF1
Other Names	CAF 1A; CAF IP60; CAF-1; CAF1; CAF1A; CAF1P60; CHAF1B; MPHOSPH7; MPP7;
Accession No.	Uniprot:Q13112
Calculated MW	Predicted band size: 61 kDa
SDS-PAGE MW	Observed band size: 70 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500-1:2000

IHC: 1:50-1:200

Images



All lanes: p60 CAF1 Rabbit mAb at 1/1k dilution

Lane 1 : Caco2 whole cell lysates
Lysates/proteins at 20 µg per lane.


Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution

Predicted band size: 61 kDa

Observed band size: 70 kDa

Exposure time: 7 seconds




All lanes:p60 CAF1 Rabbit mAb at 1/1k dilution

Lane 1 : Wild-type HeLa cell lysate

Lane 2 :p60 CAF1 Rabbit mAb knockdown HeLa cell lysate

Lysates/proteins at 20 µg per lane.



Formalin-fixed, paraffin-embedded human ovarian cancer tissue stained for p60 CAF1 using 59846 at 1/100 dilution in immunohistochemical analysis.

Background

Complex that is thought to mediate chromatin assembly in DNA replication and DNA repair. Assembles histone octamers onto replicating DNA in vitro.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.