

MST1 Rabbit mAb

Catalog No: #59848

Package Size: #59848-1 50ul #59848-2 100ul


Description

Product Name	MST1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC FC
Species Reactivity	Human Mouse Rat
Specificity	MST1 Antibody detects endogenous levels of total MST1
Immunogen Description	A synthesized peptide derived from human MST1
Other Names	HGFL; MSP; MST1; NF15S2;
Accession No.	Uniprot:P26927
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 FC 1:50

Images



Western blot analysis of MST1 expression in HepG2 cell lysate.

Background

Probably has no proteolytic activity, since crucial characteristic of serine proteases catalytic sites are not conserved.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.