

SUN1 Rabbit mAb

Catalog No: #59852

Package Size: #59852-1 50ul #59852-2 100ul


Description

Product Name	SUN1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human
Specificity	SUN1 Antibody detects endogenous levels of total SUN1
Immunogen Description	A synthesized peptide derived from human SUN1
Other Names	Sun1; UNC84A;
Accession No.	Uniprot:O94901
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of SUN1 expression in Ramos cell lysate.

Background

Component of SUN-protein-containing multivariate complexes also called LINC complexes which link the nucleoskeleton and cytoskeleton by providing versatile outer nuclear membrane attachment sites for cytoskeletal filaments.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.