

SUN2 Rabbit mAb

Catalog No: #59853

Package Size: #59853-1 50ul #59853-2 100ul


Description

Product Name	SUN2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	SUN2 Antibody detects endogenous levels of total SUN2
Immunogen Description	A synthesized peptide derived from human SUN2
Other Names	FRIGG; UNC84B;
Accession No.	Uniprot:Q9UH99
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:100~1:500 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of SUN2 expression in Jurkat cell lysate.

Background

SUN proteins form part of the LINC complex - a protein bridge that spans the nuclear envelope linking the nucleoskeleton to the actin cytoskeleton. They are located on the inner nuclear membrane side of the complex.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.