

FKBP51 Rabbit mAb

Catalog No: #59856

Package Size: #59856-1 50ul #59856-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FKBP51 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human Rat
Specificity	FKBP51 Antibody detects endogenous levels of total FKBP51
Immunogen Description	A synthesized peptide derived from human FKBP51
Other Names	AIG6; FKBP5; FKBP51; FKBP54; P54; PPIase; PPIase FKBP5; Ptg10;Rotamase;
Accession No.	Uniprot:Q13451
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50 FC 1:50

Images



Western blot analysis of FKBP51 expression in Caco-2 cell lysate.

Background

Interacts with functionally mature heterooligomeric progesterone receptor complexes along with HSP90 and TEBP.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.