

## Nectin 2 Rabbit mAb

Catalog No: #59858

Package Size: #59858-1 50ul #59858-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                                                    |
|-----------------------|----------------------------------------------------------------------------------------------------|
| Product Name          | Nectin 2 Rabbit mAb                                                                                |
| Host Species          | Rabbit                                                                                             |
| Clonality             | Monoclonal                                                                                         |
| Isotype               | Rabbit IgG                                                                                         |
| Purification          | Affinity-chromatography                                                                            |
| Applications          | WB ICC/IF IP FC                                                                                    |
| Species Reactivity    | Human Mouse Rat                                                                                    |
| Specificity           | Nectin 2 Antibody detects endogenous levels of total Nectin 2                                      |
| Immunogen Description | A synthesized peptide derived from human Nectin 2                                                  |
| Other Names           | HVEB; PRR2; CD112; PVRR2;                                                                          |
| Accession No.         | Uniprot:Q92692                                                                                     |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage               | Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.                     |

## Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 IP 1:100 FC 1:50

## Images



Western blot analysis of Nectin 2 expression in SW480 cell lysate.

## Background

Probable cell adhesion protein.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.