

CtBP1 Rabbit mAb

Catalog No: #59861

Package Size: #59861-1 50ul #59861-2 100ul


Description

Product Name	CtBP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	CtBP1 Antibody detects endogenous levels of total CtBP1
Immunogen Description	A synthesized peptide derived from human CtBP1
Other Names	BARS; CTBP; CtBP1;
Accession No.	Uniprot:Q13363
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50

Images



Western blot analysis of CtBP1 expression in Raji cell lysate.

Background

Involved in controlling the equilibrium between tubular and stacked structures in the Golgi complex. Functions in brown adipose tissue (BAT) differentiation. Corepressor targeting diverse transcription regulators such as GLIS2. Has dehydrogenase activity.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.